

## 兔抗 NPPC 多克隆抗体

中文名称：兔抗 NPPC 多克隆抗体

英文名称：Anti-NPPC rabbit polyclonal antibody

别名：INT3, NOTCH3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NPPC

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

<b>Background:</b>	This gene encodes a preproprotein that is proteolytically processed to generate multiple protein products. These products include the cardiac natriuretic peptides CNP-53, CNP-29 and CNP-22, which belong to the natriuretic family of peptides. The encoded peptides exhibit vasorelaxation activity in laboratory animals and elevated levels of CNP-22 have been observed in the plasma of chronic heart failure patients.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NPPC
<b>Immunogen:</b>	Synthetic peptide of human NPPC

<b>Full name:</b>	natriuretic peptide C
<b>Synonyms:</b>	CNP; CNP2
<b>SwissProt:</b>	P23582
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human gasrtic cancer
<b>IHC Recommend dilution:</b>	50-200

