

TRITC 标记驴抗绵羊 IgG

中文名称: TRITC 标记驴抗绵羊 IgG

英文名称: TRITC-conjugated Donkey Anti-Sheep IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标 记 物: TRITC

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody react s with whole molecule sheep IgG. It also reacts with the light cha ins of other sheep immunoglobulins. No antibody was detected a gainst non- immunoglobulin serum proteins. The antibody has be en tested by ELISA and/or solid-phase adsorbed to ensure minima I cross-reaction with chicken, guinea pig, Syrian hamster, horse, h uman, mouse, rabbit and rat serum proteins, but it may cross-rea ct with immunoglobulins from other species.
Applications	1/50-200 for most application

