

兔抗 UBN2 多克隆抗体

中文名称：兔抗 UBN2 多克隆抗体

英文名称：Anti-UBN2 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：UBN2

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Ubinuclein 2 (UBN2) is a human ortholog of the yeast protein Hpc2p. It possesses an Hpc2 (histone promoter control 2)-related domain which is highly conserved.
Applications:	ELISA, IHC
Name of antibody:	UBN2
Immunogen:	Synthetic peptide of human UBN2
Full name:	ubinnuclein 2
SwissProt:	Q6ZU65
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	20-100

