

TRITC 标记驴抗人 IgG

- 中文名称: TRITC 标记驴抗人 IgG
- 英文名称: TRITC-conjugated Donkey Anti-Human IgG
- 相关类别: 荧光标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: TRITC

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule human IgG. It also reacts with the light chain s of other human immunoglobulins. No antibody was detected ag ainst non- immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cr oss-reaction with bovine, chicken, goat, guinea pig, Syrian hamster, horse, mouse, rabbit, rat, and sheep serum proteins, but it may cr oss-react with immunoglobulins from other species.
Applications	1/50-200 for most application