

兔抗 UBOX5 多克隆抗体

中文名称：兔抗 UBOX5 多克隆抗体

英文名称：Anti-UBOX5 rabbit polyclonal antibody

别名：UIP5; RNF37; hUIP5; UBCE7IP5

抗原：UBOX5

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

相关类别：一抗

技术规格

Background:	This gene encodes a U-box domain containing protein. The encoded protein interacts with E2 enzymes and may play a role in the ubiquitination pathway. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	UBOX5
Immunogen:	Fusion protein of human UBOX5
Full name:	U-box domain containing 5
Synonyms:	UIP5; RNF37; hUIP5; UBCE7IP5

SwissProt:	O94941
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100

