

兔抗 NPPB 多克隆抗体

中文名称：兔抗 NPPB 多克隆抗体

英文名称： Anti-NPPB rabbit polyclonal antibody

别名： NPPB; BNP; Gamma-brain natriuretic peptide

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： NPPB

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	B-type natriuretic peptide (BNP) is a cardiac neurohormone specifically secreted from the ventricles in response to volume expansion and pressure overload. The actions of BNP include natriuresis, diuresis, vasorelaxation, and inhibition of renin and aldosterone secretion. It is also found to help restore the body's salt and water balance and improve heart function. BNP is co-secreted along with a biologically inactive N-terminal fragment (NT-proBNP). Circulating BNP and NT-proBNP provide important information on cardiac dysfunction, hypervolemia, and risk in patients with severe impairment of kidney function.
Applications:	WB
Name of antibody:	NPPB
Immunogen:	Fusion protein of human NPPB

Full name:	natriuretic peptide B
Synonyms:	NPPB; BNP; Gamma-brain natriuretic peptide
SwissProt:	P16860
WB Predicted band size:	15 kDa
WB Positive control:	Human plasma tissue
WB Recommended dilution:	500-2000

