

## 兔抗 RMDN2 多克隆抗体

中文名称：兔抗 RMDN2 多克隆抗体

英文名称： Anti-RMDN2 rabbit polyclonal antibody

别名： regulator of microtubule dynamics 2; RMD2; RMD4; RMD-2; FAM82A; BLOCK18; FAM82A1; PRO34163; PYST9371

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RMDN2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

The human protein FAM82A1, also known as regulator of microtubule dynamics 2 (RMD2), was one of three proteins identified as related to a family of microtubule-associated proteins in *C. elegans*. FAM82A1 contains multiple coiled-coil domains and localizes to the spindle microtubules and spindle poles during cell division. During interphase, RMD2 localizes to the cytoplasm with protein observed in the microtubule lattice and perinuclear region.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RMDN2
<b>Immunogen:</b>	Fusion protein of human RMDN2
<b>Full name:</b>	regulator of microtubule dynamics 2
<b>Synonyms:</b>	RMD2; RMD4; RMD-2; FAM82A; BLOCK18; FAM82A1; PRO34163; PYST9371
<b>SwissProt:</b>	Q96LZ7
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human tonsil
<b>IHC Recommend dilution:</b>	50-300



