

## 兔抗 RMDN3 多克隆抗体

中文名称：兔抗 RMDN3 多克隆抗体

英文名称：Anti-RMDN3 rabbit polyclonal antibody

别名：RMD3; RMD-3; FAM82C; FAM82A2; ptpip51

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RMDN3

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Involved in cellular calcium homeostasis regulation. May participate in differentiation and apoptosis of keratinocytes. Overexpression induces apoptosis. Present at high level in epidermis and seminiferous epithelium: while basal cells in the epidermis and spermatogonia show no perceptible amount, keratinocytes of suprabasal layers and differentiating first-order spermatocytes up to spermatids exhibit high expression.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	RMDN3
<b>Immunogen:</b>	Fusion protein of human RMDN3

<b>Full name:</b>	regulator of microtubule dynamics 3
<b>Synonyms :</b>	RMD3; RMD-3; FAM82C; FAM82A2; ptpip51
<b>SwissProt:</b>	Q96TC7
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human breast cancer and human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	52 kDa
<b>WB Positive control:</b>	Mouse kidney tissue
<b>WB Recommended dilution:</b>	200-1000



