

TRITC 标记驴抗山羊 IgG

中文名称: TRITC 标记驴抗山羊 IgG

英文名称: TRITC-conjugated Donkey Anti-Goat IgG

相关类别: 荧光标记二抗

储 存: 冷冻(-20℃)

标 记 物: TRITC

技术规格

Background	Based on immunoelectrophoresis and/or ELISA, the antibody rea cts with whole molecule goat IgG. It also reacts with the light c hains of other goat immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure mi nimal cross-reaction with chicken, guinea pig, Syrian hamster, ho rse, human, mouse, rabbit and rat serum proteins, but it may cr oss-react with immunoglobulins from other species.
Applications	1/50-200 for most application