

## 兔抗 UBR2 多克隆抗体

中文名称：兔抗 UBR2 多克隆抗体

英文名称：Anti-UBR2 rabbit polyclonal antibody

别名：C6orf133; bA49A4.1; dJ242G1.1; dJ392M17.3

抗原：UBR2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

相关类别：一抗

### 技术规格

<b>Background:</b>	This gene encodes an E3 ubiquitin ligase of the N-end rule proteolytic pathway that targets proteins with destabilizing N-terminal residues for polyubiquitylation and proteasome-mediated degradation. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UBR2
<b>Immunogen:</b>	Fusion protein of human UBR2
<b>Full name:</b>	ubiquitin protein ligase E3 component n-recognin 2
<b>Synonyms:</b>	C6orf133; bA49A4.1; dJ242G1.1; dJ392M17.3
<b>SwissProt:</b>	Q8IWW8
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human colorectal cancer and human cervical cancer
<b>IHC Recommend dilution:</b>	25-100



