

兔抗 RMND5B 多克隆抗体

中文名称：兔抗 RMND5B 多克隆抗体

英文名称：Anti-RMND5B rabbit polyclonal antibody

别名：required for meiotic nuclear division 5 homolog B; GID2; GID2B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RMND5B

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	RMND5B
Immunogen:	Fusion protein of human RMND5B
Full name:	required for meiotic nuclear division 5 homolog B
Synonyms:	GID2; GID2B
SwissProt:	Q96G75
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	50-300

