

兔抗 NPY2R 多克隆抗体

中文名称：兔抗 NPY2R 多克隆抗体

英文名称：Anti-NPY2R rabbit polyclonal antibody

别名：NPY2-R

相关类称：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NPY2R

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Neuropeptide Y receptor type 2 (Y2R) is a protein that in humans is encoded by the NPY2R gene. Receptor for neuropeptide Y and peptide YY. The rank order of affinity of this receptor for pancreatic polypeptides is PYY > NPY > PYY (3-36) > NPY (2-36) > [Ile-31, Gln-34] PP > [Leu-31, Pro-34] NPY > PP, [Pro-34] PYY and NPY free acid. High levels in amygdala, corpus callosum, hippocampus and subthalamic nucleus. Also detectable in caudate nucleus, hypothalamus and substantia nigra.
Applications:	ELISA, IHC
Name of antibody:	NPY2R
Immunogen:	Synthetic peptide of human NPY2R

Full name:	neuropeptide Y receptor Y2
Synonyms :	NPY2-R
SwissProt:	P49146
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human lung cancer and Human liver cancer
IHC Recommend dilution:	15-50

