

LYPLA2 抗原（重组蛋白）

中文名称： LYPLA2 抗原（重组蛋白）

英文名称： LYPLA2 Antigen (Recombinant Protein)

别名： lysophospholipase 2; APT2; APT-2; DJ886K2.4

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 1-231 amino acids of human LYPLA2

技术规格

Full name:	lysophospholipase 2
Synonyms:	APT2; APT-2; DJ886K2.4
Swissprot:	O95372
Gene Accession:	BC017034
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Lysophospholipases are enzymes that act on biological membranes to regulate the multifunctional lysophospholipids. There are alternatively spliced transcript variants described for this gene but the full length nature is not known yet.

mlbio 梅联生物
Good elisakit producers