

## 兔抗 RNASE12 多克隆抗体

- 中文名称：兔抗 RNASE12 多克隆抗体
- 英文名称：Anti-RNASE12 rabbit polyclonal antibody
- 别名：Ribonuclease-like protein 12; Rnase12; RNS12; RNASE12
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：RNASE12
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene belongs to the pancreatic ribonuclease family. It does not exhibit any ribonuclease activity. The displayed sequence is further processed into a mature form. It belongs to secreted protein.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	RNASE12
<b>Immunogen:</b>	Fusion protein of human RNASE12
<b>Full name:</b>	ribonuclease, RNase A family, 12 (non-active)
<b>Synonyms :</b>	Ribonuclease-like protein 12; Rnase12; RNS12; RNASE12
<b>SwissProt:</b>	Q5GAN4
<b>WB Predicted band size:</b>	17 kDa

<b>WB Positive control:</b>	Skin tissue
<b>WB Recommended dilution:</b>	500-2000

