

MAGEB10 抗原（重组蛋白）

中文名称：MAGEB10 抗原（重组蛋白）

英文名称：MAGEB10 Antigen (Recombinant Protein)

别名：melanoma antigen family B, 10

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 111-310 amino acids of human MAGEB10

技术规格

Full name:	melanoma antigen family B, 10
Swissprot:	Q96LZ2
Gene Accession:	BC101308
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the B subfamily of the melanoma associated antigen protein family. The encoded protein is specifically expressed in testis and tumor cells.