

## 兔抗 UBTD1 多克隆抗体

中文名称：兔抗 UBTD1 多克隆抗体

英文名称：Anti-UBTD1 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

抗 原：UBTD1

宿 主：Rabbit

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The degradation of many proteins is carried out by the ubiquitin pathway in which proteins are targeted for degradation by covalent conjugation of the poly peptide ubiquitin. This gene encodes a protein that belongs to the ubiquitin family of proteins. The encoded protein is thought to regulate E2 ubiquitin conjugating enzymes belonging to the UBE2D family.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	UBTD1
<b>Immunogen:</b>	Synthetic peptide of humanUBTD1
<b>Full name:</b>	ubiquitin domain containing 1
<b>SwissProt:</b>	Q9HAC8
<b>ELISA Recommended dilution:</b>	1000-2000

<b>IHC positive control:</b>	Human colorectal cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	26 kDa
<b>WB Positive control:</b>	A549, Jurkat, A375 and HepG2 cells
<b>WB Recommended dilution:</b>	200-1000



