

## 兔抗 RNASET2 多克隆抗体

中文名称：兔抗 RNASET2 多克隆抗体

英文名称： Anti-RNASET2 rabbit polyclonal antibody

别名： RNASE6PL; bA514O12.3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RNASET2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This ribonuclease gene is a novel member of the Rh/T2 /S-glycoprotein class of extracellular ribonucleases. It is a single copy gene that maps to 6q27, a region associated with human malignancies and chromosomal rearrangement.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RNASET2
<b>Immunogen:</b>	Fusion protein of human RNASET2
<b>Full name:</b>	ribonuclease T2
<b>Synonyms :</b>	RNASE6PL; bA514O12.3
<b>SwissProt:</b>	O00584

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	25-100

