

## 兔抗 NROB2 多克隆抗体

中文名称：兔抗 NROB2 多克隆抗体

英文名称：Anti-NROB2 rabbit polyclonal antibody

别名：nuclear receptor subfamily 0 group B member 2 SHP; SHP1

储存：冷冻（-20℃）

抗原：NROB2'

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

**Background:**

The protein encoded by this gene is an unusual orphan receptor that contains a putative ligand-binding domain but lacks a conventional DNA-binding domain. The gene product is a member of the nuclear hormone receptor family, a group of transcription factors regulated by small hydrophobic hormones, a subset of which do not have known ligands and are referred to as orphan nuclear receptors. The protein has been shown to interact with retinoid and thyroid hormone receptors, inhibiting their ligand-dependent transcriptional activation. In addition, interaction with estrogen receptors has been demonstrated, leading to inhibition of function. Studies suggest that the protein represses nuclear hormone receptor-mediated transact

	ivation via two separate steps: competition with coactivators and the direct effects of its transcriptional repressor function.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NR0B2
<b>Immunogen:</b>	Fusion protein of human NR0B2
<b>Full name:</b>	nuclear receptor subfamily 0 group B member 2
<b>Synonyms:</b>	SHP; SHP1
<b>SwissProt:</b>	Q15466
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer; Human esophagus cancer



