

MBNL3 抗原（重组蛋白）

中文名称：MBNL3 抗原（重组蛋白）

英文名称：MBNL3 Antigen (Recombinant Protein)

别名：muscleblind-like splicing regulator 3; CHCR; MBLX; MBXL; MBLX39

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 100-250 amino acids of human MBNL3

技术规格

Full name:	muscleblind-like splicing regulator 3
Synonyms:	CHCR; MBLX; MBXL; MBLX39
Swissprot:	Q9NUK0
Gene Accession:	BC042090
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the muscleblind-like family of proteins. The encoded protein may function in regulation of alternative splicing and may play a role in the pathophysiology of myotonic dystrophy. Alternatively spliced transcript variants have been described.