

兔抗 RNF11 多克隆抗体

中文名称：兔抗 RNF11 多克隆抗体

英文名称： Anti-RNF11 rabbit polyclonal antibody

别名： CGI-123; SID1669

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RNF11

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene contains a RING-H2 finger motif, which is known to be important for protein-protein interactions. The expression of this gene has been shown to be induced by mutant RET proteins (MEN2A/MEN2B). The germline mutations in RET gene are known to be responsible for the development of multiple endocrine neoplasia (MEN).
Applications:	ELISA, IHC
Name of antibody:	RNF11
Immunogen:	Fusion protein of human RNF11

Full name:	ring finger protein 11
Synonyms :	CGI-123; SID1669
SwissProt:	Q9Y3C5
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

