

## 兔抗 NR2C2AP 多克隆抗体

中文名称：兔抗 NR2C2AP 多克隆抗体

英文名称：Anti-NR2C2AP rabbit polyclonal antibody

别名：nuclear receptor 2C2 associated protein TRA16

抗原：NR2C2AP

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	May act as a repressor of NR2C2-mediated transactivation by suppressing the binding between NR2C2/TR4 and the TR4-response element in target genes.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NR2C2AP
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	nuclear receptor 2C2 associated protein
<b>Synonyms:</b>	TRA16
<b>SwissProt:</b>	Q86WQ0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

