

## 兔抗 RNF13 多克隆抗体

中文名称：兔抗 RNF13 多克隆抗体

英文名称： Anti-RNF13 rabbit polyclonal antibody

别名： ring finger protein 13; RZF

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RNF13

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. The specific function of this gene has not yet been determined. Alternatively spliced transcript variants that encode the same protein have been reported. A pseudogene, which is also located on chromosome 3, has been defined for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RNF13
<b>Immunogen:</b>	Fusion protein of human RNF13

<b>Full name:</b>	ring finger protein 13
<b>Synonyms:</b>	RZF
<b>SwissProt:</b>	O43567
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human ovarian cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	40-200



