

兔抗 RNF19A 多克隆抗体

中文名称：兔抗 RNF19A 多克隆抗体

英文名称：Anti-RNF19A rabbit polyclonal antibody

别名：BSAS; HOX1; HOX1F

相关类别：一抗

抗原：RNF19A

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the ring between ring fingers (RBR) protein family, and the encoded protein contains two RING-finger motifs and an in between RING fingers motif. This protein is an E3 ubiquitin ligase that is localized to Lewy bodies, and ubiquitylates synphilin-1, which is an interacting protein of alpha synuclein in neurons. The encoded protein may be involved in amyotrophic lateral sclerosis and Parkinson's disease. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	RNF19A

Immunogen:	Synthetic peptide corresponding to a region derived from 102-117 amino acids of human RNF19A
Full name:	ring finger protein 19A, RBR E3 ubiquitin protein ligase
Synonyms:	RNF19
SwissProt:	Q9NV58
ELISA Recommended dilution:	5000-10000
IHC Positive control:	Human liver cancer;
IHC Recommended dilution:	10-50

