

## 兔抗 NR2E1 多克隆抗体

中文名称：兔抗 NR2E1 多克隆抗体

英文名称：Anti-NR2E1 rabbit polyclonal antibody

别名：nuclear receptor subfamily 2 group E member 1; TLL; TLX; XTLL

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NR2E1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

<b>Background:</b>	The protein encoded by this gene is an orphan receptor involved in retinal development. The encoded protein also regulates adult neural stem cell proliferation and may be involved in control of aggressive behavior. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2015]
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	NR2E1
<b>Immunogen:</b>	Fusion protein of human NR2E1
<b>Full name:</b>	nuclear receptor subfamily 2 group E member 1
<b>Synonyms:</b>	TLL; TLX; XTLL

<b>SwissProt:</b>	Q9Y466
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	100-300
<b>WB Predicted band size:</b>	43 kDa
<b>WB Positive control:</b>	Human cerebrum tissue, A172 cell, Mouse brain tissue, NIH/3T3 cell lysates
<b>WB Recommended dilution:</b>	1000-5000



