

MED4 抗原（重组蛋白）

中文名称： MED4 抗原（重组蛋白）

英文名称： MED4 Antigen (Recombinant Protein)

别名： mediator complex subunit 4; ARC36; VDRIP; DRIP36; TRAP36; HSPC126

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-270 amino acids of human MED4

技术规格

Full name:	mediator complex subunit 4
Synonyms:	ARC36; VDRIP; DRIP36; TRAP36; HSPC126
Swissprot:	Q9NPJ6
Gene Accession:	BC005189
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a component of the Mediator complex. The Mediator complex interacts with DNA-binding gene-specific transcription factors to modulate transcription by RNA polymerase II. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.