

兔抗 IL6 多克隆抗体

中文名称：兔抗 IL6 多克隆抗体

英文名称：Anti-IL6 rabbit polyclonal antibody

别名：HGF; HSF; BSF2; IL-6; IFNB2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IL6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a cytokine that functions in inflammation and the maturation of B cells. In addition, the encoded protein has been shown to be an endogenous pyrogen capable of inducing fever in people with autoimmune diseases or infections. The protein is primarily produced at sites of acute and chronic inflammation, where it is secreted into the serum and induces a transcriptional inflammatory response through interleukin 6 receptor, alpha.
Applications:	ELISA, WB
Name of antibody:	IL6

Immunogen:	Synthetic peptide of human IL6
Full name:	interleukin 6
Synonyms :	HGF; HSF; BSF2; IL-6; IFNB2
SwissProt:	P05231
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	24 kDa
WB Positive control:	Hepg2 , Hela and A549 cell
WB Recommended dilution:	500-2000

