

## 兔抗 NR3C1(Ab-226) 多克隆抗体

中文名称：兔抗 NR3C1(Ab-226) 多克隆抗体

英文名称： Anti-NR3C1(Ab-226) rabbit polyclonal antibody

别名： GR; GCR; GRL; GCCR; GCRST

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： NR3C1(Ab-226)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

<b>Background:</b>	This gene encodes glucocorticoid receptor, which can function both as a transcription factor that binds to glucocorticoid response elements in the promoters of glucocorticoid responsive genes to activate their transcription, and as a regulator of other transcription factors. This receptor is typically found in the cytoplasm, but upon ligand binding, is transported into the nucleus.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	NR3C1(Ab-226)
<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human GR around the phosphorylation site of serine 226/234/246

	(L-L-S-P-L).
<b>Full name:</b>	nuclear receptor subfamily 3, group C, member 1 (glucocorticoid receptor)
<b>Synonyms :</b>	GR; GCR; GRL; GCCR; GCRST
<b>SwissProt:</b>	P04150
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	86 kDa
<b>WB Positive control:</b>	Jurkat cells treated with EGF (200ng/ml, 15mins) lysate
<b>WB Recommended dilution:</b>	500-3000



