

兔抗 NR3C1(Phospho-Ser203) 多克隆抗体

中文名称：兔抗 NR3C1(Phospho-Ser203) 多克隆抗体

英文名称： Anti-NR3C1(Phospho-Ser203) rabbit polyclonal antibody

别名： GR; GCR; GRL; GCCR; GCRST

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： NR3C1(Phospho-Ser203)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes glucocorticoid receptor, which can function both as a transcription factor that binds to glucocorticoid response elements in the promoters of glucocorticoid responsive genes to activate their transcription, and as a regulator of other transcription factors. This receptor is typically found in the cytoplasm, but upon ligand binding, is transported into the nucleus.
Applications:	WB
Name of antibody:	NR3C1(Phospho-Ser203)
Immunogen:	Peptide sequence around phosphorylation site of Serine 203 (S-G-S(p)-P-G) derived from Human GR.

Full name:	nuclear receptor subfamily 3, group C, member 1 (glucocorticoid receptor)
Synonyms :	GR; GCR; GRL; GCCR; GCRST
SwissProt:	P04150
WB Predicted band size:	86 kDa
WB Positive control:	A549 cells lysates
WB Recommended dilution:	500-1000

