

兔抗 RNF138 多克隆抗体

中文名称：兔抗 RNF138 多克隆抗体

英文名称： Anti-RNF138 rabbit polyclonal antibody

别名： ring finger protein 138; NARF; HSD-4; STRIN; hNARF

抗原： RNF138

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. Alternatively spliced transcript variants encoding distinct isoforms have been observed.
Applications:	ELISA, IHC
Name of antibody:	RNF138
Immunogen:	Synthetic peptide of human RNF138
Full name:	ring finger protein 138
Synonyms:	NARF; HSD-4; STRIN; hNARF
SwissProt:	Q8WVD3

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human liver cancer
IHC Recommend dilution:	10-50



