

## 兔抗 IL10RA(Phospho-Tyr496) 多克隆抗体

- 中文名称：兔抗 IL10RA(Phospho-Tyr496) 多克隆抗体
- 英文名称：Anti-IL10RA(Phospho-Tyr496) rabbit polyclonal antibody
- 别名：CD210; IL10R; CD210a; CDW210A; HIL-10R; IL-10R1
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：IL10RA(Phospho-Tyr496)
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a receptor for interleukin 10. This protein is structurally related to interferon receptors. It has been shown to mediate the immunosuppressive signal of interleukin 10, and thus inhibits the synthesis of proinflammatory cytokines. This receptor is reported to promote survival of progenitor myeloid cells through the insulin receptor substrate-2/PI 3-kinase/AKT pathway.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	IL10RA(Phospho-Tyr496)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of Tyrosine496 (K-G-Y(p)-L-K) derived from Human IL-10R alpha.

<b>Full name:</b>	interleukin 10 receptor, alpha
<b>Synonyms :</b>	CD210; IL10R; CD210a; CDW210A; HIL-10R; IL-10R1
<b>SwissProt:</b>	Q13651
<b>WB Predicted band size:</b>	63 kDa
<b>WB Positive control:</b>	Huvec cells lysates
<b>WB Recommended dilution:</b>	500-1000

