

兔抗 RNF148 多克隆抗体

中文名称: 兔抗 RNF148 多克隆抗体

英文名称: Anti-RNF148 rabbit polyclonal antibody

别名: ring finger protein 148

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: RNF148

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	RNF148 is a single-pass membrane protein. It contains 1 PA (protease associated) domain and 1 RING-type zinc finger. The exact function of RNF148 remains unknown.
Applications:	ELISA, IHC
Name of antibody:	RNF148
Immunogen:	Fusion protein of human RNF148
Full name:	ring finger protein 148
SwissProt:	Q8N7C7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain

IHC Recommend dilution:

50-300

