

MLLT6 抗原（重组蛋白）

中文名称：MLLT6 抗原（重组蛋白）

英文名称：MLLT6 Antigen (Recombinant Protein)

别名：MLLT6, PHD finger domain containing; AF17

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 894-1093 amino acids of human MLLT6

技术规格

Full name:	MLLT6, PHD finger domain containing
Synonyms:	AF17
Swissprot:	P55198
Gene Accession:	BC007237
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	AF17 was discovered as part of a genetic rearrangement in the MLL (mixed lineage leukemia) gene associated with acute myeloid leukemia. AF17 is also known as MLLT6.