

兔抗 NR4A1 多克隆抗体

中文名称：兔抗 NR4A1 多克隆抗体

英文名称：Anti-NR4A1 rabbit polyclonal antibody

别名：HMR; N10; TR3; NP10; GFRP1; NAK-1; NGFIB; NUR77

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NR4A1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the steroid-thyroid hormone-retinoid receptor superfamily. Expression is induced by phytohemagglutinin in human lymphocytes and by serum stimulation of arrested fibroblasts. The encoded protein acts as a nuclear transcription factor. Translocation of the protein from the nucleus to mitochondria induces apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	NR4A1
Immunogen:	Fusion protein of human NR4A1
Full name:	nuclear receptor subfamily 4, group A, member 1
Synonyms :	HMR; N10; TR3; NP10; GFRP1; NAK-1; NGFIB; NUR77

SwissProt:	P22736
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	64 kDa
WB Positive control:	PC3 cells
WB Recommended dilution:	200-1000

