

## 兔抗 IL10 多克隆抗体

- 中文名称：兔抗 IL10 多克隆抗体
- 英文名称：Anti-IL10 rabbit polyclonal antibody
- 别名：CSIF; TGIF; GVHDS; IL-10; IL10A
- 相关类别：一抗
- 储存：冷冻（-20℃）避光
- 宿主：Rabbit
- 抗原：IL10
- 反应种属：Human, Mouse, Rat, Chicken, Pig, Rabbit, Sheep, Guinea Pig
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

**Background:**

The protein encoded by this gene is a cytokine produced primarily by monocytes and to a lesser extent by lymphocytes. This cytokine has pleiotropic effects in immunoregulation and inflammation. It down-regulates the expression of Th1 cytokines, MHC class II Ags, and costimulatory molecules on macrophages. It also enhances B cell survival, proliferation, and antibody production. This cytokine can block NF-kappa B activity, and is involved in the regulation of the JAK-STAT signaling pathway. Knockout studies in mice suggested the function of this cytokine as an essential immunoregulator in the intestinal tract. Mutations in this gene are

	associated with an increased susceptibility to HIV-1 infection and rheumatoid arthritis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IL10
<b>Immunogen:</b>	KLH conjugated synthetic peptide derived from human IL-10 C-terminus
<b>Full name:</b>	interleukin 10
<b>Synonyms :</b>	CSIF; TGIF; GVHDS; IL-10; IL10A
<b>SwissProt:</b>	P22301
<b>ELISA Recommended dilution:</b>	500-1000