

兔抗 RNF207 多克隆抗体

中文名称: 兔抗 RNF207 多克隆抗体

英文名称: Anti-RNF207 rabbit polyclonal antibody

别名: ring finger protein 207; C1orf188

抗原: RNF207

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Plays a role in cardiac repolarization possibly by stabilizing membrane expression of the potassium channel KCNH2/HERG, or by assisting its synthesis, folding or export from the endoplasmic reticulum, in a heat shock protein-dependent manner.
Applications:	ELISA, IHC
Name of antibody:	RNF207
Immunogen:	Synthetic peptide of human RNF207
Full name:	ring finger protein 207
Synonyms:	C1orf188
SwissProt:	Q6ZRF8

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and Human tonsil
IHC Recommend dilution:	40-200



