

## 兔抗 UGCG 多克隆抗体

中文名称：兔抗 UGCG 多克隆抗体

英文名称： Anti-UGCG rabbit polyclonal antibody

别名： GCS; GLCT1

抗原： UGCG

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

Glycosphingolipids (GSLs) are a group of membrane components that contain lipid and sugar moieties. They are present in essentially all animal cells and are believed to have important roles in various cellular processes. UDP-glucose ceramide glucosyltransferase catalyzes the first glycosylation step in glycosphingolipid biosynthesis. The product, glucosylceramide, is the core structure of more than 300 GSLs. UGCG is widely expressed and transcription is upregulated during keratinocyte differentiation.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UGCG
<b>Immunogen:</b>	Fusion protein of human UGCG
<b>Full name:</b>	UDP-glucose ceramide glucosyltransferase
<b>Synonyms :</b>	GCS; GLCT1
<b>SwissProt:</b>	Q16739
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200

