

兔抗 NR5A1 多克隆抗体

中文名称：兔抗 NR5A1 多克隆抗体

英文名称：Anti-NR5A1 rabbit polyclonal antibody

别名：nuclear receptor subfamily 5 group A member 1; ELP; SF1; FTZ1; POF7; SF-1; AD4BP; FTZF1; SPGF8; SRXX4; SRXY3; hSF-1

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：NR5A1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a transcriptional activator involved in sex determination. The encoded protein binds DNA as a monomer. Defects in this gene are a cause of XY sex reversal with or without adrenal failure as well as adrenocortical insufficiency without ovarian defect.
Applications:	ELISA, IHC
Name of antibody:	NR5A1
Immunogen:	Fusion protein of human NR5A1
Full name:	nuclear receptor subfamily 5 group A member 1
Synonyms:	ELP; SF1; FTZ1; POF7; SF-1; AD4BP; FTZF1; SPGF8; SRXX4; SRXY3; hSF-1

SwissProt:	Q13285
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	150-300

