

兔抗 NREP 多克隆抗体

中文名称：兔抗 NREP 多克隆抗体

英文名称： Anti-NREP rabbit polyclonal antibody

别名： P311; PTZ17; SEZ17; D4S114; C5orf13; PRO1873

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NREP

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	NREP may have roles in neural function. Ectopic expression augments motility of gliomas. Promotes also axonal regeneration. May also have functions in cellular differentiation. Induces differentiation of fibroblast into myofibroblast and myofibroblast ameboid migration. Increases retinoic-acid regulation of lipid-droplet biogenesis. Down-regulates the expression of TGFB1 and TGFB2 but not of TGFB3. May play a role in the regulation of alveolar generation.
Applications:	ELISA, WB, IHC
Name of antibody:	NREP
Immunogen:	Synthetic peptide of human NREP
Full name:	neuronal regeneration related protein
Synonyms:	P311; PTZ17; SEZ17; D4S114; C5orf13; PRO1873

SwissProt:	Q16612
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human esophagus cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	8 kDa
WB Positive control:	A172 cell and Rat heart tissue lysates
WB Recommended dilution:	200-1000



