

兔抗 NRF1 多克隆抗体

中文名称：兔抗 NRF1 多克隆抗体

英文名称：Anti-NRF1 rabbit polyclonal antibody

别名：ALPHA-PAL

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NRF1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that homodimerizes and functions as a transcription factor which activates the expression of some key metabolic genes regulating cellular growth and nuclear genes required for respiration, heme biosynthesis, and mitochondrial DNA transcription and replication. The protein has also been associated with the regulation of neurite outgrowth. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB
Name of antibody:	NRF1
Immunogen:	Synthetic peptide of human NRF1
Full name:	nuclear respiratory factor 1
Synonyms :	ALPHA-PAL

SwissProt:	Q16656
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	54 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	500-2000

