

兔抗 NRG1 多克隆抗体

中文名称：兔抗 NRG1 多克隆抗体

英文名称：Anti-NRG1 rabbit polyclonal antibody

别名：GGF; HGL; HRG; NDF; ARIA; GGF2; HRG1; HRGA; SMDF; MST131; MSTP131; NRG1-IT2

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：NRG1

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a membrane glycoprotein that mediates cell-cell signaling and plays a critical role in the growth and development of multiple organ systems. An extraordinary variety of different isoforms are produced from this gene through alternative promoter usage and splicing. These isoforms are expressed in a tissue-specific manner and differ significantly in their structure, and are classified as types I, II, III, IV, V and VI. Dysregulation of this gene has been linked to diseases such as cancer, schizophrenia, and bipolar disorder (BPD).
Applications:	ELISA, IHC
Name of antibody:	NRG1

Immunogen:	Synthetic peptide of human NRG1
Full name:	neuregulin 1
Synonyms:	GGF; HGL; HRG; NDF; ARIA; GGF2; HRG1; HRGA; SMDF; MST131; MSTP131; NRG1-IT2
SwissProt:	Q02297
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and human thyroid cancer
IHC Recommend dilution:	25-100



