

## 兔抗 UGGT2 多克隆抗体

中文名称: 兔抗 UGGT2 多克隆抗体

英文名称: Anti-UGGT2 rabbit polyclonal antibody

别名: UGT2; HUGT2; UGCGL2

储存: 冷冻 (-20℃)

抗原: UGGT2

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	UDP-glucose:glycoprotein glucosyltransferase (UGT) is a soluble protein of the endoplasmic reticulum (ER) that selectively reglucosylates unfolded glycoproteins, thus providing quality control for protein transport out of the ER.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UGGT2
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	UDP-glucose glycoprotein glucosyltransferase 2
<b>Synonyms :</b>	UGT2; HUGT2; UGCGL2
<b>SwissProt:</b>	Q9NYU1

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	25-100

