

兔抗 NRG4 多克隆抗体

中文名称：兔抗 NRG4 多克隆抗体

英文名称： Anti-NRG4 rabbit polyclonal antibody

别名： HRG4

相关类别： 一抗

抗原： NRG4

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Neuregulin 4 also known as NRG4 is a member of the neuregulin protein family which in humans is encoded by the NRG4 gene. Low affinity ligand for the ERBB4 tyrosine kinase receptor. Concomitantly recruits ERBB1 and ERBB2 coreceptors, resulting in ligand-stimulated tyrosine phosphorylation and activation of the ERBB receptors. Does not bind to the ERBB1, ERBB2 and ERBB3 receptors.
Applications:	ELISA, IHC
Name of antibody:	NRG4
Immunogen:	Fusion protein of human NRG4
Full name:	neuregulin 4
Synonyms :	HRG4

SwissProt:	Q8WWG1
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human lung cancer and Human cervical cancer
IHC Recommend dilution:	50-200

