

兔抗 IL13RA2 多克隆抗体

中文名称：兔抗 IL13RA2 多克隆抗体

英文名称：Anti-IL13RA2 rabbit polyclonal antibody

别名：CT19; IL-13R; IL13BP; CD213A2

抗原：IL13RA2

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is closely related to IL13RA1, a subunit of the interleukin 13 receptor complex. This protein binds IL13 with high affinity, but lacks cytoplasmic domain, and does not appear to function as a signal mediator. It is reported to play a role in the internalization of IL13.
Applications:	ELISA, IHC
Name of antibody:	IL13RA2
Immunogen:	Fusion protein of human IL13RA2
Full name:	interleukin 13 receptor, alpha 2
Synonyms :	CT19; IL-13R; IL13BP; CD213A2
SwissProt:	Q14627

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human lymphoma and human cervical cancer
IHC Recommend dilution:	15-50

