

兔抗 NRIP1 多克隆抗体

中文名称：兔抗 NRIP1 多克隆抗体

英文名称：Anti-NRIP1 rabbit polyclonal antibody

别名：nuclear receptor interacting protein 1; CAKUT3; RIP140

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NRIP1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Nuclear receptor interacting protein 1 (NRIP1) is a nuclear protein that specifically interacts with the hormone-dependent activation domain AF2 of nuclear receptors. Also known as RIP140, this protein modulates transcriptional activity of the estrogen receptor. [provided by RefSeq, Jul 2008]
Applications:	ELISA, IHC
Name of antibody:	NRIP1
Immunogen:	Fusion protein of human NRIP1
Full name:	nuclear receptor interacting protein 1
Synonyms:	CAKUT3; RIP140
SwissProt:	P48552
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human prostate cancer and Human brain
IHC Recommend dilution:	100-200

