

兔抗 NRIP3 多克隆抗体

中文名称：兔抗 NRIP3 多克隆抗体

英文名称：Anti-NRIP3 rabbit polyclonal antibody

别名：C11orf14; NY-SAR-105

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NRIP3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	NRIP3 (nuclear receptor interacting protein 3), also known as C11orf14 or NY-SAR-105, is a 241 amino acid protein that is encoded by a gene which maps to human chromosome 11. With approximately 135 million base pairs and 1,400 genes, chromosome 11 comprises approximately 4% of human genomic DNA and is considered a gene and disease association dense chromosome. The chromosome 11 encoded Atm gene is important for regulation of cell cycle arrest and apoptosis following double strand DNA breaks.
Applications:	ELISA, IHC
Name of antibody:	NRIP3
Immunogen:	Fusion protein of human NRIP3

Full name:	nuclear receptor interacting protein 3
Synonyms :	C11orf14; NY-SAR-105
SwissProt:	Q9NQ35
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

