

## 兔抗 IL17A 多克隆抗体

- 中文名称： 兔抗 IL17A 多克隆抗体
- 英文名称： Anti-IL17A rabbit polyclonal antibody
- 别名： interleukin 17A; IL17; CTLA8; IL-17; CTLA-8; IL-17A
- 相关类别： 一抗
- 抗原： IL17A
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 反应种属： Human
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a proinflammatory cytokine produced by activated T cells. This cytokine regulates the activities of NF-kappaB and mitogen-activated protein kinases. This cytokine can stimulate the expression of IL6 and cyclooxygenase-2 (PTGS2/COX-2), as well as enhance the production of nitric oxide (NO). High levels of this cytokine are associated with several chronic inflammatory diseases including rheumatoid arthritis, psoriasis and multiple sclerosis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IL17A
<b>Immunogen:</b>	Synthetic peptide of human IL17A

<b>Full name:</b>	interleukin 17A
<b>Synonyms:</b>	IL17; CTLA8; IL-17; CTLA-8; IL-17A
<b>SwissProt:</b>	Q16552
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	30-150

