

兔抗 NRSN1 多克隆抗体

中文名称: 兔抗 NRSN1 多克隆抗体

英文名称: Anti-NRSN1 rabbit polyclonal antibody

别名: VMP; p24

抗原: NRSN1

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

相关类别: 一抗

反应种属: Human Mouse Rat

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	May play an important role in neural organelle transport, and in transduction of nerve signals or in nerve growth. May play a role in neurite extension. UniProtKB P97799
Applications:	WB
Name of antibody:	NRSN1
Immunogen:	Synthesized peptide derived from internal of human NRSN1.
Full name:	neurensin 1
Synonyms :	VMP; p24
SwissProt:	Q8IZ57
WB Predicted band size:	21 kDa
WB Positive control:	Hela cells lysate

WB Recommended dilution: 500-3000

