

## 兔抗 RPL3L 多克隆抗体

中文名称：兔抗 RPL3L 多克隆抗体

英文名称： Anti-RPL3L rabbit polyclonal antibody

储 存： 冷冻（-20℃） 避光

抗 原： RPL3L

宿 主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a protein that shares sequence similarity with ribosomal protein L3. The protein belongs to the L3P family of ribosomal proteins. Unlike the ubiquitous expression of ribosomal protein genes, this gene has a tissue-specific pattern of expression, with the highest levels of expression in skeletal muscle and heart. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome.
<b>Applications:</b>	IHC, IF
<b>Name of antibody:</b>	RPL3L
<b>Immunogen:</b>	Synthesized peptide derived from internal of human RPL3L.

<b>Full name:</b>	ribosomal protein L3-like
<b>Synonyms :</b>	0
<b>SwissProt:</b>	Q92901
<b>IHC positive control:</b>	Human heart tissue
<b>IHC Recommend dilution:</b>	50-100
<b>IF positive control:</b>	MCF-7 cells
<b>IF Recommend dilution:</b>	100-500



