

兔抗 NRTN 多克隆抗体

中文名称：兔抗 NRTN 多克隆抗体

英文名称：Anti-NRTN rabbit polyclonal antibody

别名：NTN

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NRTN

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Neurturin is a member of the TGF-beta subfamily, TRN. This gene signals through RET and a GPI-linked coreceptor, and promotes survival of neuronal populations. A neurturin mutation has been described in a family with Hirschsprung Disease.
Applications:	ELISA, IHC
Name of antibody:	NRTN
Immunogen:	Synthetic peptide of human NRTN
Full name:	neurturin
Synonyms :	NTN
SwissProt:	Q99748
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and Human breast cancer

IHC Recommend dilution: 50-200

