

兔抗 UHMK1 多克隆抗体

中文名称：兔抗 UHMK1 多克隆抗体

英文名称： Anti-UHMK1 rabbit polyclonal antibody

别名： KIS; KIST; P-CIP2

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： UHMK1

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The gene encodes a serine/threonine protein kinase that promotes cell cycle progression through G1 by phosphorylation of the cyclin-dependent kinase inhibitor 1B (p27Kip1), which causes nuclear export and degradation. The encoded protein is also thought to function in the adult nervous system and the gene has been associated with schizophrenia. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB
Name of antibody:	UHMK1
Immunogen:	Synthetic peptide of human UHMK1
Full name:	U2AF homology motif (UHM) kinase 1
Synonyms :	KIS; KIST; P-CIP2

SwissProt:	Q8TAS1
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	47 kDa
WB Positive control:	Mouse heart tissue
WB Recommended dilution:	200-1000

